

Acquisition Method Report



Agilent

Trusted Answers

Acquisition Method Info

Method Name TriPhosphateOnly_9-29-20.m
Method Path D:\MassHunter\methods\Thomas\TriPhosphateOnly_9-29-20.m
Method Description Default Method

Device List

Sampler
 Binary Pump
 Column Oven
 QQQ

MS QQQ Mass Spectrometer

Ion Source AJS ESI **Tune File** D:\MassHunter\Tune\QQQ\G6470A
 \atunes.TUNE.XML
Stop Mode By StopTime **Stop Time (min)** 40
Time Filter On **Time Filter Width (min)** 0.1

Time Segments

Index	Start Time (min)	Scan Type	Ion Mode	Div Valve	Delta EMV (+)	Delta EMV (-)	Store
1	0	MRM	ESI+Agilent Jet Stream	To MS	300	300	Yes

Time Segment 1

Scan Segments

Cpd Group	Cpd Name	ISTD?	Prec Ion	MS1 Res	Prod Ion	MS2 Res	Dwell	Frag (V)	CE (V)	Cell Acc (V)	Polarity
	UMP M+3	No	328	Unit/Enh (6490)	97	Unit/Enh (6490)	10	30	14	5	Positive
	UMP	No	325	Unit/Enh (6490)	97	Unit/Enh (6490)	10	30	14	5	Positive
	Uridine M+3	No	248.2	Unit/Enh (6490)	115.96	Unit/Enh (6490)	10	30	14	5	Positive
	Uridine	No	245.2	Unit/Enh (6490)	112.96	Unit/Enh (6490)	10	30	14	5	Positive
	ATP Standard	No	542.9	Unit/Enh (6490)	158.9	Unit/Enh (6490)	10	80	50	5	Negative
	GTP	No	522.2	Unit/Enh (6490)	159	Unit/Enh (6490)	10	30	32	5	Negative
	dGTP	No	506	Unit/Enh (6490)	159	Unit/Enh (6490)	10	30	32	5	Negative
	dATP	No	490.1	Unit/Enh (6490)	159	Unit/Enh (6490)	10	35	32	5	Negative
	UTP M+3	No	486	Unit/Enh (6490)	159	Unit/Enh (6490)	10	30	26	5	Negative
	CTP M+3	No	485.16	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	36	5	Negative
	dTTP M+3	No	484	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	35	5	Negative
	UTP	No	483	Unit/Enh (6490)	159	Unit/Enh (6490)	10	30	26	5	Negative
	CTP	No	482.16	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	36	5	Negative
	dTTP	No	481	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	35	5	Negative
	dCTP M+12 Standard	No	478	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	36	5	Negative
	dUTP M+3	No	470	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	36	5	Negative
	dCTP M+3	No	469	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	35	5	Negative
	dUTP	No	467	Unit/Enh (6490)	159	Unit/Enh (6490)	10	30	36	5	Negative

Acquisition Method Report



Agilent

Trusted Answers

Cpd Group	Cpd Name	ISTD?	Prec Ion	MS1 Res	Prod Ion	MS2 Res	Dwell	Frag (V)	CE (V)	Cell Acc (V)	Polarity
	dCTP	No	466	Unit/Enh (6490)	159	Unit/Enh (6490)	10	34	35	5	Negative
	Fumaric acid d2 Std	No	117	Unit/Enh (6490)	72.8	Unit/Enh (6490)	10	50	20	5	Negative

Source Parameters

Parameter	Value (+)	Value (-)
Gas Temp (°C)	350	350
Gas Flow (l/min)	9	9
Nebulizer (psi)	50	50
SheathGasHeater	375	375
SheathGasFlow	12	12
Capillary (V)	3500	3500
VCharging	1000	1000

Scan Parameters

Data Stg	Threshold
Centroid	0

Chromatograms

Chrom Type	Label	Offset	Y-Range
TIC	TIC	0	10000000

Instrument Curves

Actual

Name: Sampler

Module: G7129B

Auxiliary

Draw Speed	100 µL/min
Eject Speed	400 µL/min
Wait Time After Draw	1.2 s
Needle Height Offset	0.0 mm

Injection

Injection Mode	Standard injection
Injection Volume	10.00 µL

High throughput

Sample Flush-Out Factor	5.0
Injection Valve to Bypass for Delay Volume Reduction	No

Overlapped Injection

Overlapped Injection Mode	Off
---------------------------	-----

Stop Time

Stoptime Mode	As pump/No limit
---------------	------------------

Post Time

Posttime Mode	Off
---------------	-----

Name: Binary Pump

Module: G7120A

Flow	0.350 mL/min
Use Solvent Types	Yes
Stroke Mode	Synchronized
Low Pressure Limit	0.00 bar
High Pressure Limit	700.00 bar
Max. Flow Ramp Up	100.000 mL/min ²
Max. Flow Ramp Down	100.000 mL/min ²
Expected Mixer	No check

Stroke A



Automatic Stroke Calculation A

Yes

Stop Time

Stoptime Mode

Time set

Stoptime

40.00 min

Post Time

Posttime Mode

Off

Solvent Composition

	Channel	Ch. 1 Solv.	Name 1	Ch2 Solv.	Name 2	Selected	Used	Percent (%)
1	A	100.0 % Water V.02		100.0 % Water V.02		Ch. 1	Yes	15.00 %
2	B	100.0 % Acetonitrile V.02		100.0 % Acetonitrile V.02		Ch. 1	Yes	85.00 %

Timetable

	Time (min)	A (%)	B (%)	Flow (mL/min)
1	5.00 min	70.00 %	30.00 %	0.350 mL/min
2	30.00 min	98.00 %	2.00 %	0.350 mL/min
3	35.00 min	98.00 %	2.00 %	0.350 mL/min
4	40.00 min	15.00 %	85.00 %	0.350 mL/min
5	45.00 min	15.00 %	85.00 %	0.350 mL/min

Name: Column Oven

Module: G7130A

Temperature Control

Temperature Control Mode

Not Controlled

Temperature Not Ready Limit

Temperature Not Ready Limit On

Yes

Temperature Not Ready Limit Value

0.8 °C

Temperature Equilibration Time

0.0 min

Enforce column for run

Enforce column for run enabled

No

Stop Time

Stoptime Mode

As pump/injector

Post Time

Posttime Mode

Off

Timetable



	Pressure (bar)
1	700.00 bar
2	700.00 bar
3	700.00 bar
4	700.00 bar
5	700.00 bar